

Culture, Urine Quantitation - TESTING UPDATE

Effective Date: December 5, 2017

Urine collected through a straight catheter will NO longer be designated "invasive". Quantitative reporting of microorganisms from invasive and non-invasive urine specimens will follow these guidelines:

CATEGORY	INVASIVE	NON-INVASIVE
Urine Type	Autopsy	Clean Catch
	Cystoscopic	Random
	Kidney	Indwelling Catheter (Foley)
	Nephrostomy	Straight Catheter
	Suprapubic	
	Ureter	
Reporting Range (CFU/mL)	100 – 500	1,000 – 10,000
	1,000 – 5,000	10,000 – 50,000
	5,000 – 10,000	50,000 – 100,000
	>10,000	>100,000

For collection guidelines, please refer to the *Specimen Collection Manual: Urine for culture* found at: <http://beaumontlaboratory.com/test-lab-directory/lab-test-details/?testid=2109>

If you have questions, please contact Client Services (1-800-551-0488, option 5).

Laboratory Test Directory: <http://beaumontlaboratory.com/test-lab-directory>.

Date submitted: November 30, 2017
Submitted by: Bobby Boyanton, MD, Medical Director, Microbiology and Molecular Pathology

Barbara Robinson-Dunn, PhD, D(ABMM), Technical Director, Microbiology

Beaumont Laboratory
Customer Service
1-800-551-0488

28050 Grand River Ave.
Farmington Hills, MI 48336

468 Cadieux Road
Grosse Pointe, MI 48230

3601 West 13 Mile Road
Royal Oak, MI 48073

44201 Dequindre Road
Troy, MI 48085

www.beaumontlaboratory.com