

### Combined SAR-CoV-2 and Flu A B Antigens Test, Point of Care, Emergency Centers

**Effective Date: April 20, 2021**

Beaumont Laboratory will deploy SARS-CoV-2 / **Flu A B antigen** Combo Point of Care test cartridges to emergency departments in the Beaumont system. These cartridges have been validated for use in the current Sofia devices. The Combo antigen test aids in the diagnosis of COVID-19 or Influenza A/B infection in **SYMPTOMATIC** patients. Providers must order the test in Epic oneChart and answer order-entry questions required by the state.

Positive SARS-CoV-2 antigen tests in symptomatic patients should be considered **true positive**, and results do not need to be confirmed with a PCR-based test. Negative results from symptomatic patients should be treated as presumptive negative and followed up with a PCR-based test for patients to be admitted to the hospital.

Positive Flu A and/or Flu B antigen tests in symptomatic patients should be considered **presumptive** and should be confirmed with a molecular polymerase chain reaction (PCR) assay.

#### **Test Information:**

*Name:* SARS-COV-2 (COVID-19) AND INFLUENZA AB ANTIGENS, POC  
*Order Code:* LAB8483  
*CPT Code:* 87428  
*Specimen:* Nasal swab (dry)

#### **If you have questions, please contact Client Services:**

- Farmington Hills, Grosse Pointe, Royal Oak and Troy: **1-800-551-0488, option 5**
- Dearborn, Taylor, Trenton and Wayne: **1-800-245-3725, option 1**

For additional information, access the Laboratory Test Directory at:

<http://beaumontlaboratory.com/test-lab-directory>

---

**Date submitted:** April 13, 2021

**Submitted by:** Steven M. Truscott, PhD, DABCC, Technical Director, Special Chemistry, STAT Lab, Point of Care Testing  
Daniel Ortiz, PhD, D(ABMM), Associate Director, Microbiology and Molecular Pathology