

New Reference Range Dilute Russell Viper Venom Test (dRVVT)

Effective Date: April 13, 2021

In an effort to standardize reference ranges across the system, there will be a change in the reference ranges for dRVVT:

Test	New Reference Ranges	Current Reference Ranges
dRVVT Screen	<48 seconds	<44 seconds
dRVVT Confirm	<33 seconds	<35 seconds

If you have questions, please contact Client Services:

- Farmington Hills, Grosse Pointe, Royal Oak and Troy: **1-800-551-0488, option 5**
- Dearborn, Taylor, Trenton and Wayne: **1-800-245-3725, option 1**

For additional information, access the Laboratory Test Directory at:

<http://beaumontlaboratory.com/test-lab-directory>

Date submitted: March 31, 2021

Submitted by: Marc D. Smith, MD, Medical Director, Coagulation