

Beaumont Laboratory

Dearborn • Farmington Hills • Grosse Pointe • Royal Oak • Taylor • Trenton • Troy • Wayne

Urine Bilirubin Confirmation Discontinuation

Effective Date: May 16, 2023

Beaumont Laboratory will no longer perform confirmation of positive urine bilirubin results due to limited clinical utility of the test. All positive test results will be reported with the following comment:

"Positive bilirubin by dipstick. Unable to exclude color interference. Suggest clinical correlation."

Synonyms	Bilirubin in Urine/Ictotest
Instructions	
Specimen Collection Criteria	Fresh urine sample less than 2 hours old at room temperature. Fresh urine sample, less than 8 hours old refrigerated.
Physician Office/Draw Site Specimen Preparation	Fresh urine sample, less than 8 hours old refrigerated.
Specimen Preparation for Courier Transport	Fresh urine sample, less than 8 hours old refrigerated.
Rejection Criteria	Sample greater than 24 hours old, or frozen sample.
Performed	Sunday through Saturday, 24 hours per day.
Reference Range	Negative
Test Methodology	Diazotization reaction
Interpretation	The presence of bilirubin in urine is an important finding in the evaluation of liver function.
CPT Code	81003

If you have questions, please contact your Beaumont Laboratory Customer Service Department:

- Farmington Hills, Grosse Pointe, Royal Oak and Troy: 800-551-0488 or 248-551-1155, Option 5
- Dearborn, Taylor, Trenton and Wayne: 800-245-3725, Option 1

Laboratory Test Directory: <http://beaumontlaboratory.com/test-lab-directory>.

Date submitted: January 18, 2023

Submitted by: Katie (Qian) Sun, PhD, Technical Director, Automated Chemistry Core
Caitlin Schein, MD, System Medical Director