

# Laboratory Approved Specimen Collection Devices & Media

**For Corewell Health Reference Laboratory East and the following East Michigan laboratories:**

Canton, Dearborn, Farmington Hills, Grosse Pointe, Lenox, Livonia, Royal Oak, Taylor, Trenton, Troy, Wayne, West Bloomfield

Photos are for reference. Always check the container label to verify; do not go by color alone.

General test listing: not all tests are listed. Please review the test catalog for more information.\*\*

## Alinity m multi-Collect Kit



DO NOT COVER WINDOW WITH LABEL

**LAB1230885** Chlamydia, Gonococcus PCR

**LAB1231782** Chlamydia, Gonococcus, & Trichomonas PCR

**LAB1230890** STI Panel PCR\*\* (CT, NG, TV, MG)

**LAB1230886** Chlamydia PCR (CT)

**LAB1230887** Gonococcus PCR (NG)

**LAB1230889** Mycoplasma genitalium PCR (MG)

**LAB1230888** Trichomonas PCR (TV)

✓ **Sources for CT, NG, TV, or MG:** Endocervical swabs, urine, male urethral swabs, or vaginal swabs (clinician-collected).

✓ **Other sources for CT, NG, or TV:** Gynecological specimens in ThinPrep® PreservCyt solution.

✓ **Other sources for CT or NG:** Oropharyngeal swabs or rectal swabs.

## ThinPrep® Pap Container



**LAB1230097** Pap Test\*\* (Indicate HPV Reflex\*\*)

**LAB1231770** HPV High Risk Screen w/ HPV Genotype Reflex

**LAB1231782** Chlamydia, Gonococcus, & Trichomonas PCR\*

**LAB1230885** Chlamydia, Gonococcus PCR\*

**LAB1230886** Chlamydia PCR\*

**LAB1230887** Gonococcus PCR\*

**LAB1230888** Trichomonas PCR\*

✓ Aseptic Control PAP Test Lubricating Jelly and Surgilube® are the only acceptable lubricating jellies.

✓ \*Testing cannot be added onto a previous completed Pap Test Order.

## BD ESwab® Regular Swab



**LAB233**

Anaerobic Culture

**LAB240**

Fungal Culture

**LAB3417**

Group B Strep, PCR

**LAB3527**

Group B Strep Penicillin

## BD ESwab® Minitip Swab



**LAB1231047** MRSA Screen PCR


**LAB1231036** Respiratory Culture

**LAB3387** Staph aureus Screen Culture

**LAB1231048** Group A Strep PCR

**LAB2111191** Wound Culture w/gram stain

## BD® Universal Viral Transport (UVT)



**LAB9520**

Bordetella pertussis, PCR

**LAB1230607**

COVID-19 PCR

**LAB3255**

Flu A/B PCR

**LAB1231266**

Herpes Simplex Virus PCR

**LAB1231053** RSV

**LAB1231071** Varicella Zoster Virus PCR

**LAB1231483** Ureaplasma and M. hominis Culture (sendout)

## Rovers® Evalyn® Brush



**LAB1231853** Self Collect HPV vaginal swab\*\*

✓ Must be collected in a healthcare setting (i.e. physician office or health clinic). A trained healthcare staff member must immediately suspend the brush into a ThinPrep vial. See lab catalog for more information.\*\*

## Polymedco FIT Kit OC-Auto® FIT



**LAB1230599** Occult Blood, CRC Screen

**LAB1230898** Occult Blood, FIT Diagnostic

✓ Label must be filled out completely with patient legal name, date of birth and collection date and time before mailing.



## BD® Molecular Swab Kit

**LAB1231839**

Vaginitis PCR Panel



## Urinalysis and Urine Culture

**Gray Urine C&S Kit**

**Yellow UA Urine Tube**

Cup for collection only



## Para-Pak® Cary-Blair

Stool Culture  
Cryptosporidium  
Giardia



## Sterile Cup

Unpreserved Stool  
Unpreserved Urine  
>1 gram Tissue Culture  
Fluid Culture  
Sputum Culture  
Fungal Culture  
(hair skin or nail)

✓ Always use **2 patient identifiers on the specimen label** and **3 identifiers on the order**.

(Full legal name and date of birth are required)

✓ Always make sure that lids to containers, tubes and vials are on correctly and tight! Leaky specimens will be rejected.

**\*\*Review Lab Test Catalog: [corewellhealth.org/labcatalog](https://corewellhealth.org/labcatalog) | Southeast Michigan**

ThinPrep® and Aptima® are registered trademarks of Hologic, Inc. BD® and ESwab® is a registered trademark of Becton, Dickinson and Company. UTM® is a registered trademark of Copan Diagnostics. OC-Auto® FIT is a registered trademark of Polymedco, LLC. Surgilube® is a registered trademark of HR Pharmaceuticals, Inc. Evalyn® Brush is a registered trademark of Rovers® Medical Devices.

March 2025 (East)