

LABORATORY

Beaumont Laboratory's tube chart shows the most commonly ordered tests with their associated specimen tube. Please confirm proper preparation and most updated tube type in our Laboratory Test Directory or call our client service department at **800-245-3725**.

SERUM SEPARATOR TUBE (SST)

Acute Hep Pl.	Bilirubin	Ferritin	HCV Ab	Iron/TIBC	Phosphorus	PSA	TSH
AST	BUN	Folic Acid	hCG Beta Quant	LD	Potassium	Renal Func. Pl.	Uric Acid
ALT	Creatinine	FSH	Hepatic Func. Pl.	Lipase	Progesterone	Rheum. Factor	Vit B-12
Amylase	CRP-HS	HBcAB, Total	HIV 1/2 Ab	Lipid Pl.	Protein Electro	RPR	Vit D, 25 Hydrox
Alk. Phos.	Comp. Met. Pl.	HBsAg	Homocysteine	LH	Protein, Total	T4, Free	
Basic Met Pl.	Electrolytes Pl.	HBsAb	HSV 1 & 2, IgG	Magnesium	PTH *	Testosterone	

LAVENDAR TUBE K₂ EDTA

BNP	Hgb A1C	Evaluation	Path Consult
CBC/Plt	Hemoglobin and	HIV Quant	Retic Count
CBC/Diff & Plt	Hematocrit	Lead Level,	Sed Rate
Factor V Leiden	Hemoglobinopathy	Whole Blood	Vit B-6 **

LT. BLUE 3.2% Sodium Citrate

Anti-Thrombin III*	Hexagonal Phase	Protein S Activity *
D-Dimer *	Phospholipid *	PTINR
DRVVT *	Inhibitor Screen *	PTT
Factor V Activity *	Protein C Activity *	

RED TUBE Plain

Beta Hydroxybutyrate	Phenobarbital	Theophylline
Carbamazepine	Phenytoin	Tobramycin
Digoxin	Phenytoin, Free	Valproic Acid
Ionized Ca	T3, Free	Vancomycin

PINK TUBE EDTA

ABO & RH Type
Ab Screen, RBC
ABO
RH Type

GREY TUBE Potassium Oxalate/ Sodium Fluoride

Glucose Tolerance Test	Glucose Tolerance 2hr, Pregnancy
Glucose Pregnancy Screen, 1hr	Glucose Tolerance 3hr, Pregnancy

GREEN TUBE Sodium Herapin

Ammonia (transport on ice)
Lactic Acid (transport on ice)

LT. GREEN TUBE Lithium Herapin Plasma Gel. Sep.

Homocysteine
Troponin I

ROYAL BLUE TUBE EDTA

Heavy Metals Screen	Mercury
Magnesium, RBC	Zinc RBC

ROYAL BLUE TUBE Plain

Aluminum	Nickel
Copper	Zinc Level
Manganese	

• frozen pour over - special preparation * protect from light