

Specimen Collection Guide

ESwab*

Regular Tip ESwab (White Cap)	
Soft Test Code	Test Name
AGSSA	Strep Screen, Group A
CXANA	Culture, Anaerobic
CXEAR	Culture, Ear
CXEYE	Culture, Eye
CXGC	Culture, GC (Female)
CXGEN	Culture, Genital
CXMDR	Culture, MDR-XDR
CXMRS	Culture, MRSA
CXRES	Culture, Respiratory (Throat, Pharynx)
CXSA	Culture, <i>Staph. aureus</i>
CXTIS	Culture, Tissue (< 1 gram)
CXVRE	Culture, VRE
CXWND	Culture, Wound Deep
CXWNS	Culture, Wound Superficial
SSA	Group A Strep by NAA without Susceptibility
SSAS	Group A Strep by NAA with Susceptibility
SSB	Group B Strep Screen
STNOR	Stain, Nocardia
WP	Wet Prep Exam



Mini Tip ESwab (Green Cap)	
Soft Test Code	Test Name
CXGC	Culture, GC (Male)
CXRES	Culture, Respiratory (Nasopharynx)
IBPNA	<i>Bordetella pertussis</i> by nucleic acid amplification



- * Do not submit the following using ESwab:
- Nucleic acid amplification tests – exception, *Bordetella pertussis* (IBPNA)
 - Sterile body fluids
 - Tissue (>1 gram)

For more information or questions, please contact Customer Service at 800-551-0488.

Beaumont

BEAUMONT LABORATORY