

Total Hemoglobin and Oxyhemoglobin in Whole Blood, Point of Care, Cath Lab

Effective Date: 10/17/2022

Beaumont Laboratory will replace the AVOX/OPL Point of Care tests in the Dearborn Cath LAB with the Radiometer ABL80 instruments for the tests below.

The AVOX/OPL Point of Care instruments will be retired.

All Total Hemoglobin and Oxyhemoglobin in Whole Blood will be tested on the ABL80s.

The ABL80s are hard wire connected and will be interfaced to automatically send acceptable results to the patient charts in EPIC.

The current and new test panels are displayed in the chart below:

AVOX/OPL	ABL80
Hemoglobin	Hemoglobin
Oxyhemoglobin	Oxyhemoglobin

Test Information:

Name: Co-Oximetry, Arterial, POC
Order Code: LAB7151
CPT Codes: 85018 (hemoglobin), 82820 (Oxyhemoglobin)
Specimen: 1.5 mL heparinized whole blood in a plastic syringe or *safePICO* syringe

Name: Co-Oximetry, Cardiac Panel, POC
Order Code: LAB7150
CPT Codes: 85018 (hemoglobin), 82820 (Oxyhemoglobin)
Specimen: 1.5 mL heparinized whole blood in a plastic syringe or *safePICO* syringe

Name: Co-Oximetry, Mixed Venous, POC
Order Code: LAB7153
CPT Codes: 85018 (hemoglobin), 82820 (Oxyhemoglobin)
Specimen: 1.5 mL heparinized whole blood in a plastic syringe or *safePICO* syringe

Name: Co-Oximetry, Venous, POC
Order Code: LAB7152
CPT Codes: 85018 (hemoglobin), 82820 (Oxyhemoglobin)
Specimen: 1.5 mL heparinized whole blood in a plastic syringe or *safePICO* syringe

Date submitted: 10/5/2022

Submitted by: Nga Yeung Tang, PhD, Technical Director, STAT Lab, Point of Care Testing
Jeremy Powers, MD, Chief Dearborn Pathology